

# UROLOGY

This catalog remains exclusive property of the company H. + H. Maslanka Chirurgische Instrumente GmbH. Reprinting, even partially, is not allowed. Technical alterations and printing errors reserved.



Autoclavable at 134°C / 273°F - 5 minutes



Sterilised with gamma rays



Sterilised with ethylene oxid



Singe use only - The instruments must not be resterilised respectively reprocessed.



Radio-opaque - The instruments are visible with gamma rays.



With port for contrast medium.



With cleaning port



With high frequency connector - Please consider the generator and instrument settings.




## Urology

	<b>page</b>
Biopsy forceps . . . . .	4-11
Foreign Body Retriever . . . . .	12-33
Scissors . . . . .	34-41
Stone basket . . . . .	42-53
Injection needle . . . . .	54-58
High frequency snare . . . . .	59-63
Cleaning brushes . . . . .	64-68
Cytology brush . . . . .	69-71
Coagulation probe . . . . .	72-75
Accessory for high-frequency . . . . .	76-78

## Urology

### Biopsy forceps



	<b>page</b>
Biopsy forceps - reusable - Product Group U020 .....	5-11



## Urology

### Biopsy forceps

#### Biopsy forceps - reusable - Product Group U020

	<b>page</b>
Biopsy forceps with oval cups - flexible . . . . .	6
Biopsy forceps with oval cups - rigid . . . . .	7
Biopsy forceps with toothed cups flexible - interlocked . . . . .	8
Biopsy forceps with toothed cups rigid - interlocked . . . . .	9
Biopsy forceps with toothed cups flexible - successive toothed . . . . .	10
Biopsy forceps with toothed cups rigid - successive toothed . . . . .	11

## Biopsy forceps with oval cups - flexible



- Oval cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
0,8 mm - 2,4 Fr.	on request	<b>S0208-040</b>	<b>S0208-060</b>
1,0 mm - 3 Fr.	<b>S0210-030</b>	<b>S0210-040</b>	<b>S0210-060</b>
1,2 mm - 3,5 Fr.	on request	<b>S0212-040</b>	<b>S0212-060</b>
1,6 mm - 5 Fr.	<b>S0316-030</b>	<b>S0316-040</b>	<b>S0316-060</b>
1,8 mm - 5,5 Fr.	on request	<b>S0318-040</b>	<b>S0318-060</b>
2,1 mm - 6 Fr.	on request	<b>S0221-040</b>	<b>S0221-060</b>
2,3 mm - 7 Fr.	on request	<b>S0222-040</b>	<b>S0222-060</b>

Required accessories: Handle and cable

Please find the appropriate article numbers on our accessory pages "accessory for high-frequency"

## Biopsy forceps with oval cups - rigid



- Oval cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
0,8 mm - 2,4 Fr.	on request	<b>RS0208-040</b>	<b>RS0208-060</b>
1,0 mm - 3 Fr.	<b>RS0210-030</b>	<b>RS0210-040</b>	<b>RS0210-060</b>
1,2 mm - 3,5 Fr.	on request	<b>RS0212-040</b>	<b>RS0212-060</b>
1,6 mm - 5 Fr.	<b>RS0316-030</b>	<b>RS0316-040</b>	<b>RS0316-060</b>
1,8 mm - 5,5 Fr.	on request	<b>RS0318-040</b>	<b>RS0318-060</b>
2,1 mm - 6 Fr.	on request	<b>RS0221-040</b>	<b>RS0221-060</b>
2,3 mm - 7 Fr.	on request	<b>RS0222-040</b>	<b>RS0222-060</b>

Required accessories: Handle and cable

Please find the appropriate article numbers on our accessory pages "accessory for high-frequency"



## Biopsy forceps with toothed cups flexible - interlocked

- Toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
1,0 mm - 3 Fr.	<b>S0410-030</b>	<b>S0410-040</b>	<b>S0410-060</b>
1,6 mm - 5 Fr.	<b>S0416-030</b>	<b>S0416-040</b>	<b>S0416-060</b>
1,8 mm - 5,5 Fr.	<b>S0418-030</b>	<b>S0418-040</b>	<b>S0418-060</b>
2,1 mm - 6,3 Fr.	<b>S0421-030</b>	<b>S0421-040</b>	<b>S0421-060</b>
2,3 mm - 7 Fr.	<b>S0422-030</b>	<b>S0422-040</b>	<b>S0422-060</b>

Required accessories: Handle and cable

Please find the appropriate article numbers on our accessory pages "accessory for high-frequency"





## Biopsy forceps with toothed cups rigid - interlocked

- Toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
1,0 mm - 3 Fr.	<b>RS0410-030</b>	<b>RS0410-040</b>	<b>RS0410-060</b>
1,6 mm - 5 Fr.	<b>RS0416-030</b>	<b>RS0416-040</b>	<b>RS0416-060</b>
1,8 mm - 5,5 Fr.	<b>RS0418-030</b>	<b>RS0418-040</b>	<b>RS0418-060</b>
2,1 mm - 6,3 Fr.	<b>RS0421-030</b>	<b>RS0421-040</b>	<b>RS0421-060</b>
2,3 mm - 7 Fr.	<b>RS0422-030</b>	<b>RS0422-040</b>	<b>RS0422-060</b>

Required accessories: Handle and cable  
Please find the appropriate article numbers on our accessory pages "accessory for high-frequency"



## Biopsy forceps with toothed cups flexible - successive toothed

- Toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
1,6 mm - 5 Fr.	<b>S04416-030</b>	<b>S04416-040</b>	<b>S04416-060</b>
2,1 mm - 6,3 Fr.	<b>S04421-030</b>	<b>S04421-040</b>	<b>S04421-060</b>
2,3 mm - 7 Fr.	<b>S04422-030</b>	<b>S04422-040</b>	<b>S04422-060</b>

Required accessories: Handle and cable

Please find the appropriate article numbers on our accessory pages "accessory for high-frequency"



## Biopsy forceps with toothed cups rigid - successive toothed

- Toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
1,6 mm - 5 Fr.	<b>RS04416-030</b>	<b>RS04416-040</b>	<b>RS04416-060</b>
2,1 mm - 6,3 Fr.	<b>RS04421-030</b>	<b>RS04421-040</b>	<b>RS04421-060</b>
2,3 mm - 7 Fr.	<b>RS04422-030</b>	<b>RS04422-040</b>	<b>RS04422-060</b>

Required accessories: Handle and cable  
Please find the appropriate article numbers on our accessory pages "accessory for high-frequency"

## Urology

### Foreign Body Retriever



**page**

Foreign Body Retriever - reusable - Product Group U022 ..... 13-33




## Urology

### Foreign Body Retriever

#### Foreign Body Retriever - reusable - Product Group U022

	<b>page</b>
Foreign Body Retriever serrated jaws - 4,8 mm flexible - successive toothed . . . . .	14
Foreign Body Retriever serrated jaws - 4,8 mm rigid - successive toothed . . . . .	15
Foreign Body Retriever serrated jaws 6,6 mm long - flexible - interlocked . . . . .	16
Foreign Body Retriever serrated jaws 6,6 mm - rigid - interlocked toothed . . . . .	17
Foreign Body Retriever - flexible - interlocked toothed . . . . .	18
Foreign Body Retriever - rigid - interlocked toothed . . . . .	19
Foreign Body Retriever - one sided flexible - interlocked toothed . . . . .	20
Foreign Body Retriever - one sided rigid - interlocked toothed . . . . .	21
Foreign Body Retriever - flexible . . . . .	22
Foreign Body Retriever - rigid . . . . .	23
Foreign Body Retriever for stone-street - flexible . . . . .	24
Foreign Body Retriever for stone-street - rigid . . . . .	25
DJ - Forceps - flexible . . . . .	26
DJ - Forceps - rigid . . . . .	27
Rat tooth forceps - flexible . . . . .	28
Rat tooth forceps - rigid . . . . .	29
Foreign Body Retriever - flexible . . . . .	30
Foreign Body Retriever - rigid . . . . .	31
Forceps for stones and fragments - rigid . . . . .	32
Forceps for stones and fragments - type curved . . . . .	33




## Foreign Body Retriever serrated jaws - 4,8 mm flexible - successive toothed

- Cross toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
1,2 mm - 3,5 Fr.	on request	<b>S04612-040</b>	<b>S04612-060</b>
1,6 mm - 5 Fr.	on request	<b>S04616-040</b>	<b>S04616-060</b>



## Foreign Body Retriever serrated jaws - 4,8 mm rigid - successive toothed

- Cross toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
1,2 mm - 3,5 Fr.	on request	<b>RS04612-040</b>	<b>RS04612-060</b>
1,6 mm - 5 Fr.	on request	<b>RS04616-040</b>	<b>RS04616-060</b>



## Foreign Body Retriever serrated jaws 6,6 mm long - flexible - interlocked

- Cross toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
1,6 mm - 5 Fr.	on request	<b>R04316-040</b>	<b>R04316-060</b>
2,1 mm - 6 Fr.	on request	<b>R04321-040</b>	<b>R04321-060</b>






## Foreign Body Retriever serrated jaws 6,6 mm - rigid - interlocked toothed

- Cross toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
1,6 mm - 5 Fr.	on request	<b>RS04316-040</b>	<b>RS04316-060</b>
2,1 mm - 6 Fr.	on request	<b>RS04321-040</b>	<b>RS04321-060</b>




## Foreign Body Retriever - flexible - interlocked toothed

- Cross toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
1,2 mm - 3,5 Fr.	on request	<b>S05212-040</b>	<b>S05212-060</b>
2,3 mm - 7 Fr.	on request	<b>S05222-040</b>	<b>S05222-060</b>



## Foreign Body Retriever - rigid - interlocked toothed

- Cross toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
1,2 mm - 3,5 Fr.	on request	<b>RS05212-040</b>	<b>RS05212-060</b>
2,3 mm - 7 Fr.	on request	<b>RS05222-040</b>	<b>RS05222-060</b>



## Foreign Body Retriever - one sided flexible - interlocked toothed

- Cross toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length
	40 cm
2,3 mm - 7 Fr.	<b>S05322-040</b>




## Foreign Body Retriever - one sided rigid - interlocked toothed

- Cross toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length
	40 cm
2,3 mm	RS05322-040



## Foreign Body Retriever - flexible


- Cross toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
1,0 mm - 3 Fr.	<b>S1810-030</b>	<b>S1810-040</b>	<b>S1810-060</b>
1,2 mm - 3,5 Fr.	on request	<b>S1812-040</b>	<b>S1812-060</b>
1,6 mm - 5 Fr.	<b>S1816-030</b>	<b>S1816-040</b>	<b>S1816-060</b>
1,8 mm - 5,5 Fr.	on request	<b>S1818-040</b>	<b>S1818-060</b>
2,4 mm - 7 Fr.	<b>S1822-030</b>	<b>S1822-040</b>	<b>S1822-060</b>

Required accessories: Handle and cable

Please find the appropriate article numbers on our accessory pages "accessory for high-frequency"



## Foreign Body Retriever - rigid

- Cross toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
1,0 mm - 3 Fr.	<b>RS1810-030</b>	<b>RS1810-040</b>	<b>RS1810-060</b>
1,2 mm - 3,5 Fr.	on request	<b>RS1812-040</b>	<b>RS1812-060</b>
1,6 mm - 5 Fr.	<b>RS1816-030</b>	<b>RS1816-040</b>	<b>RS1816-060</b>
1,8 mm - 5,5 Fr.	on request	<b>RS1818-040</b>	<b>RS1818-060</b>
2,4 mm - 7 Fr.	<b>RS1822-030</b>	<b>RS1822-040</b>	<b>RS1822-060</b>

Required accessories: Handle and cable  
Please find the appropriate article numbers on our accessory pages "accessory for high-frequency"



## Foreign Body Retriever for stone-street - flexible

- Long cross toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length
1,6 mm - 5 Fr.	60 cm <b>S04116-060</b>






## Foreign Body Retriever for stone-street - rigid

- Long cross toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length
1,6 mm - 5 Fr.	60 cm RS04116-060




## DJ - Forceps - flexible

- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length	
	40 cm	60 cm
2,3 mm - 7 Fr.	<b>S08322-040</b>	<b>S08322-060</b>
3,0 mm - 9 Fr.	<b>S08330-040</b>	<b>S08330-060</b>



## DJ - Forceps - rigid

- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length	
	40 cm	60 cm
2,3 mm - 7 Fr.	<b>RS08322-040</b>	<b>RS08322-060</b>
3,0 mm - 9 Fr.	<b>RS08330-040</b>	<b>RS08330-060</b>



## Rat tooth forceps - flexible

- 2:1 teeth
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length	
	40 cm	60 cm
1,6 mm - 5 Fr.	<b>S08416-040</b>	<b>S08416-060</b>
2,3 mm - 7 Fr.	<b>S08422-040</b>	<b>S08422-060</b>



## Rat tooth forceps - rigid

- 2:1 teeth
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length	
	40 cm	60 cm
1,6 mm - 5 Fr.	<b>RS08416-040</b>	<b>RS08416-060</b>
2,3 mm - 7 Fr.	<b>RS08422-040</b>	<b>RS08422-060</b>

## Foreign Body Retriever - flexible



- Forked jaws with 1:2 tooth
- Jaws not totally retractable
- Handle with cleaning port



Instrument diameter	Working length
	60 cm
1,6 mm - 5 Fr.	<b>S1016-060</b>
1,8 mm - 5,5 Fr.	<b>S1018-060</b>

## Foreign Body Retriever - rigid



- Forked jaws with 1:2 tooth
- Jaws not totally retractable
- Handle with cleaning port



Instrument diameter	Working length
	60 cm
1,6 mm - 5 Fr.	<b>RS1016-060</b>
1,8 mm - 5,5 Fr.	<b>RS1018-060</b>

## Forceps for stones and fragments - rigid



- Fenestrated
- Ground spiral
- handle with cleaning port



Instrument diameter	Working length
3,0 mm - 9 Fr.	40 cm RS05130-040





## Forceps for stones and fragments - type curved

- Type curved - rigid
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length
3,0 mm - 9 Fr.	30 cm <b>RS18130-030</b>

## Urology

### Scissors



	<b>page</b>
Scissors - reusable - Product Group U024 .....	35-41




## Urology

### Scissors

#### Scissors - reusable - Product Group U024

	<b>page</b>
Scissors one-sided - blunt tips flexible . . . . .	36
Scissors one-sided - blunt tips rigid . . . . .	37
Scissors double-sided - blunt tips rigid . . . . .	38
Scissors one-sided - pointed tips flexible . . . . .	39
Scissors one-sided - pointed rigid . . . . .	40
Hook Scissors, single action . . . . .	41




## Scissors one-sided - blunt tips flexible

- One-sided movable scissors
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
1,0 mm - 3 Fr.	<b>S09103-030</b>	<b>S09103-040</b>	<b>S09103-060</b>
1,6 mm - 5 Fr.	<b>S09163-030</b>	<b>S09163-040</b>	<b>S09163-060</b>
2,1 mm - 6 Fr.	<b>S09213-030</b>	<b>S09213-040</b>	<b>S09213-060</b>
2,3 mm - 7 Fr.	<b>S09223-030</b>	<b>S09223-040</b>	<b>S09223-060</b>



## Scissors one-sided - blunt tips rigid

- One-sided movable scissors
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
1,0 mm - 3 Fr.	<b>RS09103-030</b>	<b>RS09103-040</b>	<b>RS09103-060</b>
1,6 mm - 5 Fr.	<b>RS09163-030</b>	<b>RS09163-040</b>	<b>RS09163-060</b>
2,1 mm - 6 Fr.	<b>RS09213-030</b>	<b>RS09213-040</b>	<b>RS09213-060</b>
2,3 mm - 7 Fr.	<b>RS09223-030</b>	<b>RS09223-040</b>	<b>RS09223-060</b>



## Scissors double-sided - blunt tips rigid

- Double-sided movable scissors
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length	
	30 cm	40 cm
1,6 mm - 5 Fr.	<b>RS0916-030</b>	<b>RS0916-040</b>



## Scissors one-sided - pointed tips flexible

- One-sided movable scissors
- Ground spiral
- Handle with cleaning port
- Pointed tips



Instrument diameter	Working length	
	30 cm	40 cm
1,6 mm - 5 Fr.	<b>S09162-030</b>	<b>S09162-040</b>
2,3 mm - 7 Fr.	<b>S09222-030</b>	<b>S09222-040</b>



## Scissors one-sided - pointed rigid

- One-sided movable scissors
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length	
	30 cm	40 cm
1,6 mm - 5 Fr.	<b>RS09162-030</b>	<b>RS09162-040</b>
2,3 mm - 7 Fr.	<b>RS09222-030</b>	<b>RS09222-040</b>



## Hook Scissors, single action



- One-sided movable scissors
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length
2,3 mm - 7 Fr.	70 cm RS09224

## Urology

### Stone basket



	<b>page</b>
Stone basket - single use - Product Group U031 .....	43-47
Stone basket - reusable - Product Group U030 .....	48-53



## Urology

### Stone basket

#### Stone basket - single use - Product Group U031

	<b>page</b>
Stone basket 4-fold turned - single use Ø 1,0 mm - 3 Fr. - Nitinol, form stable . . . . .	44
Stone basket 4-fold turned - single use Ø 1,3 mm - 4 Fr. - Nitinol, form stable . . . . .	45
Stone basket 4-fold turned - single use Ø 1,65 mm - 5 Fr. . . . .	46
Stone basket 4-fold turned - single use Ø 1,8 mm - 5,5 Fr. . . . .	47



### Stone basket 4-fold turned - single use

**Ø 1,0 mm - 3 Fr. - Nitinol, form stable**

For working channel Ø 1,2 mm and higher

- Twisted construction
- 4 wires - Nitinol, form stable
- Fixed handle



Basket type	Basket diameter	Working length	
		70 cm	90 cm
Nitinol wire	15 mm	<b>E154.410-070</b>	<b>E154.410-090</b>



### Stone basket 4-fold turned - single use

Ø 1,3 mm - 4 Fr. - Nitinol, form stable

For working channel Ø 1,5 mm and higher

- Twisted construction
- 4 wires - Nitinol, form stable
- Fixed handle



Basket type	Basket diameter	Working length	
		70 cm	90 cm
Nitinol wire	15 mm	E154.413-070	E154.413-090



## Stone basket 4-fold turned - single use

Ø 1,65 mm - 5 Fr.

For working channel Ø 1,8 mm and higher

- Twisted construction
- 4 wires
- With Luer Lok port for contrast medium and cleaning
- Fixed handle



Basket type	Basket diameter	Working length
		60 cm
Monofile wire	15 mm	E150164-060
	20 mm	E150165-060

Also available in nitinol! Please enquire.



## Stone basket 4-fold turned - single use

**Ø 1,8 mm - 5,5 Fr.**

For working channel Ø 2,0 mm and higher

- Twisted construction
- 4 wires
- With Luer Lok port for contrast medium and cleaning
- Fixed handle



Basket type	Basket diameter	Working length
		60 cm
Monofile wire	15 mm	E150184-060
	20 mm	E150185-060

Also available in nitinol! Please enquire.



## Urology

### Stone basket

#### Stone basket - reusable - Product Group U030

	<b>page</b>
Stone basket 4-fold turned Ø 1,0 mm - 3 Fr. ....	49
Stone basket 4-fold turned Ø 1,65 mm - 5 Fr. ....	50
Stone basket 4-fold turned Ø 1,8 mm - 5,5 Fr. ....	51
Stone basket 6-fold turned Ø 1,65 mm - 5 Fr. ....	52
Stone basket 6-fold turned Ø 1,8 mm - 5,5 Fr. ....	53





### Stone basket 4-fold turned

**Ø 1,0 mm - 3 Fr.**

For working channel Ø 1,2 mm and higher

- Twisted construction
- 4 wires
- With Luer Lok port for contrast medium and cleaning
- Can be taken apart for cleaning
- With handle



Basket Type	Basket diameter	Working length	
		40 cm	60 cm
Monofile wire	15 mm	<b>150104-040</b>	<b>150104-060</b>



### Stone basket 4-fold turned

Ø 1,65 mm - 5 Fr.

For working channel Ø 1,8 mm and higher

- Twisted construction
- 4 wires
- With Luer Lok port for contrast medium and cleaning
- Can be taken apart for cleaning
- With handle



Basket type	Basket diameter	Working length	
		40 cm	60 cm
Monofile wire	15 mm	150164-040	150164-060
	20 mm	150165-040	150165-060
Multifile wire	15 mm	160164-040	160164-060
	23 mm	160165-040	160165-060



## Stone basket 4-fold turned

**Ø 1,8 mm - 5,5 Fr.**

For working channel Ø 2,0 mm and higher

- Twisted construction
- 4 wires
- With Luer Lok port for contrast medium and cleaning
- Can be taken apart for cleaning
- With handle



Basket type	Basket diameter	Working length	
		40 cm	60 cm
Monofile wire	15 mm	<b>150184-040</b>	<b>150184-060</b>
	20 mm	<b>150185-040</b>	<b>150185-060</b>
Multifile wire	15 mm	<b>160184-040</b>	<b>160184-060</b>
	23 mm	<b>160185-040</b>	<b>160185-060</b>



### Stone basket 6-fold turned

**Ø 1,65 mm - 5 Fr.**

For working channel Ø 1,8 mm and higher

- Twisted construction
- 6 wires
- With Luer Lok port for contrast medium and cleaning
- Can be taken apart for cleaning
- With handle



Basket type	Basket diameter	Working length	
		40 cm	60 cm
Monofile wire	15 mm	<b>152164-040</b>	<b>152164-060</b>
	20 mm	<b>152165-040</b>	<b>152165-060</b>
Multifile wire	15 mm	<b>162164-040</b>	<b>162164-060</b>
	23 mm	<b>162165-040</b>	<b>162165-060</b>



### Stone basket 6-fold turned

Ø 1,8 mm - 5,5 Fr.

For working channel Ø 2,0 mm and higher


- Twisted construction
- 6 wires
- With Luer Lok port for contrast medium and cleaning
- Can be taken apart for cleaning
- With handle



Basket type	Basket diameter	Working length	
		40 cm	60 cm
Monofile wire	15 mm	<b>152184-040</b>	<b>152184-060</b>
	20 mm	<b>152185-040</b>	<b>152185-060</b>
Multifile wire	15 mm	<b>162184-040</b>	<b>162184-060</b>
	23 mm	<b>162185-040</b>	<b>162185-060</b>

## Urology

### Injection needle



	<b>page</b>
Injection needle - single use - Product Group U001 .....	55-56
Injection needle - reusable - Product Group U006 .....	57-58



## Urology

### Injection needle

#### Injection needle - single use - Product Group U001

	<b>page</b>
Injection needle with needle - 5,5 Fr. single use .....	56



## Injection needle with needle - 5,5 Fr. single use

For working channels Ø 2,0 mm and higher

- Needle diameter 0,7 mm
- Needle length 4 mm or 6 mm
- Sterile packed



Needle length	Working Length	Type
		Plastic tube
4 mm	60 cm	<b>ET2518.74-060</b>
6 mm	60 cm	<b>ET2518.76-060</b>





## Urology

### Injection needle

#### Injection needle - reusable - Product Group U006

	page
Injection needle with needle - 5,5 Fr. ....	58



## Injection needle with needle - 5,5 Fr.


- Needle diameter 0,7 mm
- 60 cm working length



Needle length	Model	Type
		Plastic tube
10 mm	Needle retractable	2518-50/10
10 mm	Fixed	2518-50/10F

## Urology

### High frequency snare



**page**

High frequency snare - single use - Product Group U005 ..... 60-63



## Urology

### High frequency snare

#### High frequency snare - single use - Product Group U005

	<b>page</b>
Symmetrical High Frequency Snare Ø 1,65 mm - 5 Fr. - single use . . . . .	61
Symmetrical High Frequency Snare Ø 1,8 mm - 5,5Fr. - single use . . . . .	62
Symmetrical High Frequency Snare Ø 2,35 mm - 7 Fr. - single use . . . . .	63



## Symmetrical High Frequency Snare Ø 1,65 mm - 5 Fr. - single use

For working channels Ø 1,8 mm and higher


- Symmetrical snare-opening
- Available in two loop sizes
- With Handle



HF connection type	Working length	Loop size	
		20 mm	25 mm
4 mm socket	60 cm	<b>E341615-060</b>	<b>E341616-060</b>
Olympus	60 cm	<b>E341615A-060</b>	<b>E341616A-060</b>

Required accessories: Cable

Please find the appropriate article numbers on our accessory pages "accessory for high-frequency"



## Symmetrical High Frequency Snare Ø 1,8 mm - 5,5Fr. - single use

For working channels Ø 2,0 mm and higher

- Symmetrical snare-opening
- Available in two loop sizes
- With Handle



HF connection type	Working length	Loop size	
		20 mm	25 mm
4 mm socket	60 cm	<b>E341815-060</b>	<b>E341816-060</b>
Olympus	60 cm	<b>E341815A-060</b>	<b>E341816A-060</b>

Required accessories: Cable

Please find the appropriate article numbers on our accessory pages "accessory for high-frequency"



## Symmetrical High Frequency Snare Ø 2,35 mm - 7 Fr. - single use

For working channels Ø 2,45 mm and higher

- Symmetrical snare-opening
- Available in two loop sizes
- With Handle



HF connection type	Working length	Loop size	
		20 mm	25 mm
4 mm socket	60 cm	<b>E342215-060</b>	<b>E342216-060</b>
Olympus	60 cm	<b>E342215A-060</b>	<b>E342216A-060</b>

Required accessories: Cable

Please find the appropriate article numbers on our accessory pages "accessory for high-frequency"

## Urology

### Cleaning brushes



	<b>page</b>
Cleaning brush single use - Product Group U004 .....	65-66
Cleaning brush - reusable - Product Group U004 .....	67-68





## Urology

### Cleaning brushes

#### Cleaning brush single use - Product Group U004

page

Cleaning brush - single use . . . . . 66

## Cleaning brush - single use



- double-sided brush
- plastic tip for protection
- different working channels



	Working length	For working channel			
		Ø 1,0 mm	Ø 1,5 mm	Ø 2,0 mm	Ø 2,8 mm
Double-sided	60 cm	<b>EB4110-060</b>	<b>EB4115-060</b>	<b>EB4118-060</b>	<b>EB4122-060</b>



## Urology

### Cleaning brushes

#### Cleaning brush - reusable - Product Group U004

page

Cleaning brush .....68

## Cleaning brush




- With single oder double-sided brush
- Flexible distal tip
- For different working channels



	Working length	For working channel			
		Ø 1,2 mm and higher	Ø 1,65 mm and higher	Ø 2,0 mm and higher	Ø 2,6 mm and higher
Single-sided	40 cm	<b>41101-040</b>	<b>41161-040</b>	<b>41181-040</b>	<b>41221-040</b>
	60 cm	<b>41101-060</b>	<b>41161-060</b>	<b>41181-060</b>	<b>41221-060</b>
Double-sided	40 cm	<b>BS41101-040</b>	<b>BS41161-040</b>	<b>BS41181-040</b>	<b>BS41221-040</b>
	60 cm	<b>BS41101-060</b>	<b>BS41161-060</b>	<b>BS41181-060</b>	<b>BS41221-060</b>

## Urology

### Cytology brush



**page**

Cytology brush for single use - Product Group U003 . . . . . 70-71



## Urology

### Cytology brush

#### Cytology brush for single use - Product Group U003

page

Cytology brush for single-use .....71



## Cytology brush for single-use

- With plastic tube
- Brush totally retractable
- With ball for tube sealing
- Single sterile packed
- Packing unit 10 pieces



Working length	For working channel
	Ø 1,2 mm and higher
120 cm	<b>E4210-120</b>

## Urology

### Coagulation probe



**page**

Coagulation probe -reusable - Product Group U007 .....73-75





## Urology

### Coagulation probe

#### Coagulation probe -reusable - Product Group U007

	<b>page</b>
Coagulation probe Ø 1,65 mm - 5 Fr. ....	74
Coagulation probe Ø 1,8 mm - 5,5 Fr. ....	75



## Coagulation probe Ø 1,65 mm - 5 Fr.

For working channels Ø 1,8 mm and higher

- With or without Axial flushing port
- Flexible



Spiral type	Working length	Art. No.
With flushing port	40 cm	<b>8016-040</b>
	60 cm	<b>8016-060</b>
Without flushing port	40 cm	<b>80161-040</b>
	60 cm	<b>80161-060</b>

Required accessories: Cable

Please find the appropriate article numbers on our accessory pages "accessory for high-frequency"



## Coagulation probe Ø 1,8 mm - 5,5 Fr.

For working channels Ø 2,0 mm and higher

- With or without Axial flushing port
- Flexible



Spiral type	Working length	Art.No.
With flushing port	40 cm	<b>8018-040</b>
	60 cm	<b>8018-060</b>
Without flushing port	40 cm	<b>80181-040</b>
	60 cm	<b>80181-060</b>

Required accessories: Cable

Please find the appropriate article numbers on our accessory pages "accessory for high-frequency"

## Urology

### Accessory for high-frequency



page

High frequency cable - Product Group 0033 ..... 77-78



## Urology

### Accessory for high-frequency

#### High frequency cable - Product Group 0033

	page
High frequency cable .....	78

## High frequency cable



- 4 mm connector (Martin, Erbe, Storz)
- 5 mm connector (Erbe)
- 6 mm connector (Olympus)
- 8 mm connector (Vallylab)



Generator side	Handle side	
	4 mm socket	Olympus
Martin Erbe (ACC) Storz	<b>6400</b>	<b>6400A</b>
Erbe (ICC)	<b>6404</b>	<b>6404A</b>
Olympus	<b>6402</b>	<b>6402A</b>
Vallylab	<b>6407</b>	<b>6407A</b>