







Precise removal

Reusable polypectomy & mucosectomy snares

MASLANKA reusable polypectomy & mucosectomy snares

Benefit from the versatility of our snares




| Symbol | Meaning |
|---|-----------------------|
|  | Without insulated tip |
|  | With insulated tip |
|  | Multifilament |
|  | Monofilament |
|  | Rotatable |

Flexibility & functionality in focus.


Our reusable snares, offer highest user comfort, whether with **multifilament or monofilament cutting wire**. The **non-insulated tip** enables selective coagulation of the removal site by the user after removal of polyps, without having to change instruments.

You can also benefit from our reusable **multifilament** snares **with insulated tip**. This prevents perforation during extension of the snare or coagulation during ablation of the polyp under poor visibility.



Without insulated tip

| Article number | Working channel | Snare opening | Working length | PU |
|----------------|-----------------|---------------|----------------|---------|
| 34181/55-230 | Ø 2,0 mm | - 20 mm | 230 cm | 1 piece |
| 34181/65-230 | Ø 2,0 mm | - 25 mm | 230 cm | 1 piece |
| 34181/75-230 | Ø 2,0 mm | - 35 mm | 230 cm | 1 piece |
| 34221/55-230 | Ø 2,8 mm | - 20 mm | 230 cm | 1 piece |
| 34221/65-230 | Ø 2,8 mm | - 25 mm | 230 cm | 1 piece |
| 34221/75-230 | Ø 2,8 mm | - 35 mm | 230 cm | 1 piece |




Without insulated tip, with fixed handle (Olympus HF connector)

| Article number | Working channel | Snare opening | Working length | PU |
|----------------|-----------------|---------------|----------------|---------|
| K34221/55A-230 | Ø 2,8 mm | - 20 mm | 230 cm | 1 piece |
| K34221/65A-230 | Ø 2,8 mm | - 25 mm | 230 cm | 1 piece |
| K34221/75A-230 | Ø 2,8 mm | - 35 mm | 230 cm | 1 piece |


Rotatable, without insulated tip

| Article number | Working channel | Snare opening | Working length | PU |
|----------------|-----------------|---------------|----------------|---------|
| D34221/55-230 | Ø 2,8 mm | - 20 mm | 230 cm | 1 piece |
| D34221/65-230 | Ø 2,8 mm | - 25 mm | 230 cm | 1 piece |
| D34221/75-230 | Ø 2,8 mm | - 35 mm | 230 cm | 1 piece |



With insulated tip

| Article number | Working channel | Snare opening | Working length | PU |
|----------------|-----------------|---------------|----------------|---------|
| 34180/55-230 | Ø 2,0 mm | - 20 mm | 230 cm | 1 piece |
| 34180/65-230 | Ø 2,0 mm | - 25 mm | 230 cm | 1 piece |
| 34180/75-230 | Ø 2,0 mm | - 35 mm | 230 cm | 1 piece |
| 34220/55-230 | Ø 2,8 mm | - 20 mm | 230 cm | 1 piece |
| 34220/65-230 | Ø 2,8 mm | - 25 mm | 230 cm | 1 piece |
| 34220/75-230 | Ø 2,8 mm | - 35 mm | 230 cm | 1 piece |





With insulated tip, with fixed handle (Olympus HF connector)

| Article number | Working channel | Snare opening | Working length | PU |
|----------------|-----------------|---------------|----------------|---------|
| K34220/55A-230 | Ø 2,8 mm | - 20 mm | 230 cm | 1 piece |
| K34220/65A-230 | Ø 2,8 mm | - 25 mm | 230 cm | 1 piece |
| K34220/75A-230 | Ø 2,8 mm | - 35 mm | 230 cm | 1 piece |

Rotatable, with insulated tip


| Article number | Working channel | Snare opening | Working length | PU |
|----------------|-----------------|---------------|----------------|---------|
| D34220/55-230 | Ø 2,8 mm | - 20 mm | 230 cm | 1 piece |
| D34220/65-230 | Ø 2,8 mm | - 25 mm | 230 cm | 1 piece |
| D34220/75-230 | Ø 2,8 mm | - 35 mm | 230 cm | 1 piece |





Reusable
Symmetrical monofilament polypectomy snare

Without insulated tip

| Article number | Working channel | Snare opening | Working length | PU |
|----------------|-----------------|---------------|----------------|---------|
| 34221M/35-230 | Ø 2,8 mm | - 10 mm | 230 cm | 1 piece |
| 34221M/55-230 | Ø 2,8 mm | - 20 mm | 230 cm | 1 piece |
| 34221M/65-230 | Ø 2,8 mm | - 25 mm | 230 cm | 1 piece |
| 34221M/75-230 | Ø 2,8 mm | - 30 mm | 230 cm | 1 piece |





Symmetrical multifilament mucosectomy snare

Available only with fixed handle.



Without insulated tip (Olympus HF connector)


| Article number | Working channel | Snare opening | Working length | PU |
|----------------|-----------------|---------------|----------------|---------|
| K34226/75A-230 | Ø 2,8 mm | ~ 30 mm | 230 cm | 1 piece |



Asymmetrical multifilament polypectomy snare

Without insulated tip

| Article number | Working channel | Snare opening | Working length | PU |
|----------------|-----------------|---------------|----------------|---------|
| 3218-180 | Ø 2,0 mm | - 15 mm | 180 cm | 1 piece |
| 3218-230 | Ø 2,0 mm | - 15 mm | 230 cm | 1 piece |
| 3322-180 | Ø 2,8 mm | - 25 mm | 180 cm | 1 piece |
| 3322-230 | Ø 2,8 mm | - 25 mm | 230 cm | 1 piece |



With insulated tip

| Article number | Working channel | Snare opening | Working length | PU |
|----------------|-----------------|---------------|----------------|---------|
| 3022-230 | Ø 2,8 mm | - 15 mm | 230 cm | 1 piece |
| 3122-230 | Ø 2,8 mm | - 25 mm | 230 cm | 1 piece |



Without insulated tip, with fixed handle (Olympus HF connector)


| Article number | Working channel | Snare opening | Working length | PU |
|----------------|-----------------|---------------|----------------|---------|
| K3218A-180 | Ø 2,0 mm | - 15 mm | 180 cm | 1 piece |
| K3218A-230 | Ø 2,0 mm | - 15 mm | 230 cm | 1 piece |
| K3322A-180 | Ø 2,8 mm | - 25 mm | 180 cm | 1 piece |
| K3322A-230 | Ø 2,8 mm | - 25 mm | 230 cm | 1 piece |

With insulated tip, with fixed handle (Olympus HF connector)

| Article number | Working channel | Snare opening | Working length | PU |
|----------------|-----------------|---------------|----------------|---------|
| K3022A-230 | Ø 2,8 mm | - 15 mm | 230 cm | 1 piece |
| K3122A-230 | Ø 2,8 mm | - 25 mm | 230 cm | 1 piece |


Rotatable, without isolated tip


| Article number | Working channel | Snare opening | Working length | PU |
|----------------|-----------------|---------------|----------------|---------|
| D3322-230 | Ø 2,8 mm | - 25 mm | 230 cm | 1 piece |



Rotatable with insulated tip

| Article number | Working channel | Snare opening | Working length | PU |
|----------------|-----------------|---------------|----------------|---------|
| D3122-230 | Ø 2,8 mm | - 25 mm | 230 cm | 1 piece |





134°C Autoclave 273°F Reusable

Asymmetrical multifilament polypectomy snare




Asymmetrical multifilament polypectomy snare with push handle Reusable



With push handle, without isolated tip (Olympus HF connector)

| Article number | Working channel | Snare opening | Working length | PU |
|----------------|-----------------|---------------|----------------|---------|
| 6522.40A-180 | Ø 2,8 mm | - 40 mm | 180 cm | 1 piece |
| 6522.40A-230 | Ø 2,8 mm | - 40 mm | 230 cm | 1 piece |
| 6522.60A-180 | Ø 2,8 mm | - 60 mm | 180 cm | 1 piece |
| 6522.60A-230 | Ø 2,8 mm | - 60 mm | 230 cm | 1 piece |


Scale 1:1




With push handle, with isolated tip (Olympus HF connector)

| Article number | Working channel | Snare opening | Working length | PU |
|----------------|-----------------|---------------|----------------|---------|
| 6522.401A-180 | Ø 2,8 mm | - 40 mm | 180 cm | 1 piece |
| 6522.401A-230 | Ø 2,8 mm | - 40 mm | 230 cm | 1 piece |
| 6522.601A-180 | Ø 2,8 mm | - 60 mm | 180 cm | 1 piece |
| 6522.601A-230 | Ø 2,8 mm | - 60 mm | 230 cm | 1 piece |

Scale 1:1






 **Reusable**
Hexagonal multifilament polypectomy snare



Without insulated tip

| Article number | Working channel | Snare opening | Working length | PU |
|----------------|-----------------|---------------|----------------|---------|
| 37180/65-230 | Ø 2,0 mm | - 25 mm | 230 cm | 1 piece |
| 37220/55-230 | Ø 2,8 mm | - 20 mm | 230 cm | 1 piece |
| 37220/65-230 | Ø 2,8 mm | - 25 mm | 230 cm | 1 piece |
| 37220/75-230 | Ø 2,8 mm | - 35 mm | 230 cm | 1 piece |




Without insulated tip, with fixed handle (Olympus HF connector)

| Article number | Working channel | Snare opening | Working length | PU |
|----------------|-----------------|---------------|----------------|---------|
| K37180/65A-230 | Ø 2,0 mm | - 25 mm | 230 cm | 1 piece |
| K37220/55A-230 | Ø 2,8 mm | - 20 mm | 230 cm | 1 piece |
| K37220/65A-230 | Ø 2,8 mm | - 25 mm | 230 cm | 1 piece |
| K37220/75A-230 | Ø 2,8 mm | - 35 mm | 230 cm | 1 piece |




Rotatable, without isolated tip

| Article number | Working channel | Snare opening | Working length | PU |
|----------------|-----------------|---------------|----------------|---------|
| D37180/65-230 | Ø 2,0 mm | - 25 mm | 230 cm | 1 piece |
| D37220/55-230 | Ø 2,8 mm | - 20 mm | 230 cm | 1 piece |
| D37220/65-230 | Ø 2,8 mm | - 25 mm | 230 cm | 1 piece |
| D37220/75-230 | Ø 2,8 mm | - 35 mm | 230 cm | 1 piece |



With insulated tip

| Article number | Working channel | Snare opening | Working length | PU |
|----------------|-----------------|---------------|----------------|---------|
| 38180/65-230 | Ø 2,0 mm | - 25 mm | 230 cm | 1 piece |
| 38220/55-230 | Ø 2,8 mm | - 20 mm | 230 cm | 1 piece |
| 38220/65-230 | Ø 2,8 mm | - 25 mm | 230 cm | 1 piece |
| 38220/75-230 | Ø 2,8 mm | - 35 mm | 230 cm | 1 piece |



Luer lock connector for quick and easy mounting

Simple mounting to the cleaning machine




Reusable

Tube attachment for non-rotatable models

With insulated tip, with fixed handle (Olympus HF connector)

| Article number | Working channel | Snare opening | Working length | PU |
|----------------|-----------------|---------------|----------------|---------|
| K38180/65A-230 | Ø 2,0 mm | - 25 mm | 230 cm | 1 piece |
| K38220/55A-230 | Ø 2,8 mm | - 20 mm | 230 cm | 1 piece |
| K38220/65A-230 | Ø 2,8 mm | - 25 mm | 230 cm | 1 piece |
| K38220/75A-230 | Ø 2,8 mm | - 35 mm | 230 cm | 1 piece |



Luer lock connector for quick and easy mounting of the handle

Rotatable

Simple mounting to the cleaning machine




Reusable

Tube attachment for rotatable models

Rotatable, with insulated tip

| Article number | Working channel | Snare opening | Working length | PU |
|----------------|-----------------|---------------|----------------|---------|
| D38180/65-230 | Ø 2,0 mm | - 25 mm | 230 cm | 1 piece |
| D38220/55-230 | Ø 2,8 mm | - 20 mm | 230 cm | 1 piece |
| D38220/65-230 | Ø 2,8 mm | - 25 mm | 230 cm | 1 piece |
| D38220/75-230 | Ø 2,8 mm | - 35 mm | 230 cm | 1 piece |



Reusable Handles




134°C Autoclave 273°F Reusable

Handle with Olympus HF connector for tube attachment


Suitable for standard handle & handle for rotatable versions.

| Article number | Connection | PU |
|----------------|-------------------------------------|---------|
| 6203LL | Olympus HF connector for 4 mm cable | 1 piece |



134°C Autoclave 273°F Reusable
Fixed handle

More products from Maslanka



Our product portfolio includes products for gastroenterology, cardiology and urology, as well as specialist solutions for surgery.

Our product catalogue contains a detailed breakdown of our products. Please do not hesitate to contact us if you require any further information.

Gastroenterology
Cardiology
Urology
Neurology
Customer-specific developments

Reusable Polypectomy & mucosectomy snares



H. + H. Maslanka
Chirurgische Instrumente GmbH
Stockacherstraße 172
78532 Tuttlingen / Germany

Tel.: +49 7461 9607 - 0
Fax: +49 7461 - 47 32
info@maslanka.de
www.maslanka.de

06/2022 | EN

CE 0483