



Even Application

Reusable laparoscopic applicator
for use with Protein-Thrombin-Matrix

MASLANKA Reusable laparoscopic applicator

No metal abrasion due to the use of high-quality, stainless steels

Safe and even application of pressure thanks to the ergonomic plastic handle

Can't find the product you want?
Talk to us to find out more.



Consists of two components:
A non-reflective stainless steel trocar and a stainless steel mandrel with plastic handle for better cleaning

Slight overhang for complete application



Reusable

Laparoscopic applicator
for use with Protein-Thrombin-Matrix
including cleaning brush

Article number	Working channel from diam.	Working length	PU
285040-35	5,0 mm	35 cm	1