



Disposable



# Complete protection

Disposable cytology brush

MASLANKA Disposable cytology brush

Closure of the working channel with the ball head **prevents cross-contamination**

Ball head as **protection against damage** in the endoscope

Stable shaft for **optimised force transfer**

**Short distance** between the head and brush for direct contact

Ergonomic handle

**Swivel brush** thanks to rotating ring

## Disposable cytology brush

Our disposable product.

Article number	Working channel	Brush diameter	Brush length	Working length	PU
E4218-120	Ø 2,0 mm	Ø 2,2 mm	15 mm	120cm	10 pcs
E4218-180	Ø 2,0 mm	Ø 2,2 mm	15 mm	180 cm	10 pcs
E4222-230	Ø 2,8 mm	Ø 2,2 mm	15 mm	230 cm	10 pcs