



Injection needles  
for urology




H. + H. Maslanka  
Chirurgische Instrumente GmbH  
Stockacherstraße 172  
78532 Tuttlingen / Germany

Tel.: +49 (0) 7461 9607 - 0  
Fax: +49 (0) 7461 - 47 32  
info@maslanka.de  
www.maslanka.de

# Spot on

Injection needles for urology


CE 0483 03/2020 | EN



134°C Autoclave 273°F **Multi-use**  
 Injection needles for urology, **retractable**

Our multi-use product.

Article number	Needle diameter	Needle length	Working channel	Working length	PU	
2518-50/10	Ø 0.7 mm	10 mm	2.0 mm	60 cm	10 pcs	needle retractable



134°C Autoclave 273°F **Multi-use**  
 Injection needles for urology, **fixed**

Our multi-use product.

Article number	Needle diameter	Needle length	Working channel	Working length	PU	
2518-50/10F	Ø 0.7 mm	10 mm	2.0 mm	60 cm	10 pcs	needle fixed