






Flexibility as a strength

Reusable papillotomes

MASLANKA Reusable papillotomes



Reusable
Papillotome, **single-lumen**

With fixed handle


Article number	Working channel from	Design	Cutting length	Working length	SU
K4101825	Diam. 2,0 mm	Filiform	~ 25 mm	200 cm	1 unit
K4121825		Filiform, with precut	~ 25 mm		1 unit



Reusable
Papillotome for 0,035 inch **guide wire, single lumen**

With fixed handle

Article number	Working channel from	Design	Cutting length	Working length	SU
K4162230	Ø 2,8 mm	Filiform, distally tapered to diam. 1.85 mm	~ 30 mm	200 cm	1 unit
K4172230		Filiform, cutting at 11 o'clock	~ 30 mm		1 unit




134°C
273°F
Reusable

Papillotome for 0,035 inch **guide wire, double-lumen**

With fixed handle

Article number	Working channel from	Design	Cutting length	Working length	SU
K5102220	Ø 2,8 mm	Filiform, distally tapered to diam. 2.1 mm	~ 20 mm	200 cm	1 unit
K5102225		Filiform, distally tapered to diam. 2.1 mm	~ 25 mm		1 unit
K5102230		Filiform, distally tapered to diam. 2.1 mm	~ 30 mm		1 unit
K5012225		Standard, with precut, distally tapered to diam. 2.1 mm	~ 25 mm		1 unit




134°C
273°F
Reusable

Straight Papillotome **single-lumen**

With fixed handle

Type	Article number	For a working channel from	Design	Working length	SU
Needle type	K46018	Diam. 2,0 mm	6 mm	200 cm	1 unit
Blade type	K47018		6 mm		1 unit



Diam 2,8 mm





Reusable Papillotome for **B-II-stomachs**

With fixed handle

Article number	For a working channel from	Cutting length	Working length	SU
K49022	Ø 3,2 mm	~ 20 mm	200 cm	1 Stück

More products from Maslanka



Products for gastroenterology, cardiology and urology as well as special solutions for surgery are all part of our product range.

You will find a detailed overview of our products in our product catalogue. Please do not hesitate to contact us if you have any further questions.

Gastroenterology

Cardiology

Urology

Neurology

Customer-specific developments

Reusable Papillotomes



H. + H. Maslanka
Chirurgische Instrumente GmbH
Stockacherstraße 172
78532 Tuttlingen / Germany

Tel.: +49 7461 9607 - 0
Fax: +49 7461 - 47 32
info@maslanka.de
www.maslanka.de

09/2020 | DE

CE 0483