


# Precise removal

Reusable polypectomy & mucosectomy snares

MASLANKA reusable polypectomy & mucosectomy snares

# Benefit from the versatility of our snares




| Symbol | Meaning               |
|--------|-----------------------|
|        | Without insulated tip |
|        | With insulated tip    |
|        | Multifilament         |
|        | Monofilament          |
|        | Rotatable             |

## Flexibility & functionality in focus.

Our reusable snares, offer highest user comfort, whether with **multifilament or monofilament cutting wire**. The **non-insulated tip** enables selective coagulation of the removal site by the user after removal of polyps, without having to change instruments.

You can also benefit from our reusable **multifilament** snares **with insulated tip**. This prevents perforation during extension of the snare or coagulation during ablation of the polyp under poor visibility.



134°C  
Autoclave  
273°F

Reusable

Symmetrical multifilament polypectomy snare

Without insulated tip

| Article number | Working channel | Snare opening | Working length | PU      |
|----------------|-----------------|---------------|----------------|---------|
| 34181/55-230   | Ø 2,0 mm        | - 20 mm       | 230 cm         | 1 piece |
| 34181/65-230   | Ø 2,0 mm        | - 25 mm       | 230 cm         | 1 piece |
| 34181/75-230   | Ø 2,0 mm        | - 35 mm       | 230 cm         | 1 piece |
| 34221/55-230   | Ø 2,8 mm        | - 20 mm       | 230 cm         | 1 piece |
| 34221/65-230   | Ø 2,8 mm        | - 25 mm       | 230 cm         | 1 piece |
| 34221/75-230   | Ø 2,8 mm        | - 35 mm       | 230 cm         | 1 piece |

Scale 1:1

Without insulated tip, with fixed handle (Olympus HF connector)


| Article number | Working channel | Snare opening | Working length | PU      |
|----------------|-----------------|---------------|----------------|---------|
| K34221/55A-230 | Ø 2,8 mm        | - 20 mm       | 230 cm         | 1 piece |
| K34221/65A-230 | Ø 2,8 mm        | - 25 mm       | 230 cm         | 1 piece |
| K34221/75A-230 | Ø 2,8 mm        | - 35 mm       | 230 cm         | 1 piece |

Rotatable, without insulated tip


| Article number | Working channel | Snare opening | Working length | PU      |
|----------------|-----------------|---------------|----------------|---------|
| D34221/55-230  | Ø 2,8 mm        | - 20 mm       | 230 cm         | 1 piece |
| D34221/65-230  | Ø 2,8 mm        | - 25 mm       | 230 cm         | 1 piece |
| D34221/75-230  | Ø 2,8 mm        | - 35 mm       | 230 cm         | 1 piece |

Scale 1:1


| With insulated tip |                 |               |                |         |
|--------------------|-----------------|---------------|----------------|---------|
| Article number     | Working channel | Snare opening | Working length | PU      |
| 34180/55-230       | Ø 2,0 mm        | - 20 mm       | 230 cm         | 1 piece |
| 34180/65-230       | Ø 2,0 mm        | - 25 mm       | 230 cm         | 1 piece |
| 34180/75-230       | Ø 2,0 mm        | - 35 mm       | 230 cm         | 1 piece |
| 34220/55-230       | Ø 2,8 mm        | - 20 mm       | 230 cm         | 1 piece |
| 34220/65-230       | Ø 2,8 mm        | - 25 mm       | 230 cm         | 1 piece |
| 34220/75-230       | Ø 2,8 mm        | - 35 mm       | 230 cm         | 1 piece |



| With insulated tip, with fixed handle (Olympus HF connector) |                 |               |                |         |
|--|-----------------|---------------|----------------|---------|
| Article number   | Working channel | Snare opening | Working length | PU      |
| K34220/55A-230   | Ø 2,8 mm        | - 20 mm       | 230 cm         | 1 piece |
| K34220/65A-230   | Ø 2,8 mm        | - 25 mm       | 230 cm         | 1 piece |
| K34220/75A-230   | Ø 2,8 mm        | - 35 mm       | 230 cm         | 1 piece |




| Rotatable, with insulated tip |                 |               |                |         |
|-------------------------------|-----------------|---------------|----------------|---------|
| Article number                | Working channel | Snare opening | Working length | PU      |
| D34220/55-230                 | Ø 2,8 mm        | - 20 mm       | 230 cm         | 1 piece |
| D34220/65-230                 | Ø 2,8 mm        | - 25 mm       | 230 cm         | 1 piece |
| D34220/75-230                 | Ø 2,8 mm        | - 35 mm       | 230 cm         | 1 piece |



Reusable

## Symmetrical monofilament polypectomy snare

| Without insulated tip |                 |               |                |         |
|-----------------------|-----------------|---------------|----------------|---------|
| Article number        | Working channel | Snare opening | Working length | PU      |
| 34221M/35-230         | Ø 2,8 mm        | - 10 mm       | 230 cm         | 1 piece |
| 34221M/55-230         | Ø 2,8 mm        | - 20 mm       | 230 cm         | 1 piece |
| 34221M/65-230         | Ø 2,8 mm        | - 25 mm       | 230 cm         | 1 piece |
| 34221M/75-230         | Ø 2,8 mm        | - 30 mm       | 230 cm         | 1 piece |





Reusable



## Symmetrical multifilament mucosectomy snare

Available only with fixed handle.



### Without insulated tip (Olympus HF connector)

| Article number | Working channel | Snare opening | Working length | PU      |
|----------------|-----------------|---------------|----------------|---------|
| K34226A/75-230 | Ø 2,8 mm        | ~ 30 mm       | 230 cm         | 1 piece |



Reusable

## Asymmetrical multifilament polypectomy snare

Without insulated tip

| Article number | Working channel | Snare opening | Working length | PU      |
|----------------|-----------------|---------------|----------------|---------|
| 3218-180       | Ø 2,0 mm        | ~ 15 mm       | 160 cm         | 1 piece |
| 3218-230       | Ø 2,0 mm        | ~ 15 mm       | 230 cm         | 1 piece |
| 3322-180       | Ø 2,8 mm        | ~ 25 mm       | 160 cm         | 1 piece |
| 3322-230       | Ø 2,8 mm        | ~ 25 mm       | 230 cm         | 1 piece |

Scale 1:1



With insulated tip

| Article number | Working channel | Snare opening | Working length | PU      |
|----------------|-----------------|---------------|----------------|---------|
| 3022-230       | Ø 2,8 mm        | ~ 15 mm       | 230 cm         | 1 piece |
| 3122-230       | Ø 2,8 mm        | ~ 25 mm       | 230 cm         | 1 piece |

Scale 1:1



Without insulated tip, with fixed handle (Olympus HF connector)

| Article number | Working channel | Snare opening | Working length | PU      |
|----------------|-----------------|---------------|----------------|---------|
| K3218A-180     | Ø 2,0 mm        | ~ 15 mm       | 160 cm         | 1 piece |
| K3218A-230     | Ø 2,0 mm        | ~ 15 mm       | 230 cm         | 1 piece |
| K3322A-180     | Ø 2,8 mm        | ~ 25 mm       | 160 cm         | 1 piece |
| K3322A-230     | Ø 2,8 mm        | ~ 25 mm       | 230 cm         | 1 piece |

With insulated tip, with fixed handle (Olympus HF connector)

| Article number | Working channel | Snare opening | Working length | PU      |
|----------------|-----------------|---------------|----------------|---------|
| K3022A-230     | Ø 2,8 mm        | ~ 15 mm       | 230 cm         | 1 piece |
| K3122A-230     | Ø 2,8 mm        | ~ 25 mm       | 230 cm         | 1 piece |

Rotatable, without isolated tip

| Article number | Working channel | Snare opening | Working length | PU      |
|----------------|-----------------|---------------|----------------|---------|
| D3322-230      | Ø 2,8 mm        | ~ 25 mm       | 230 cm         | 1 piece |

Scale 1:1





Rotatable with insulated tip

| Article number | Working channel | Snare opening | Working length | PU      |
|----------------|-----------------|---------------|----------------|---------|
| D3122-230      | Ø 2,8 mm        | ~ 25 mm       | 230 cm         | 1 piece |

Scale 1:1








#### With push handle, without isolated tip (Olympus HF connector)

| Article number | Working channel | Snare opening | Working length | PU      |
|----------------|-----------------|---------------|----------------|---------|
| 6522.40A-180   | Ø 2,8 mm        | ~ 40 mm       | 180 cm         | 1 piece |
| 6522.40A-230   | Ø 2,8 mm        | ~ 40 mm       | 230 cm         | 1 piece |
| 6522.60A-180   | Ø 2,8 mm        | ~ 60 mm       | 180 cm         | 1 piece |
| 6522.60A-230   | Ø 2,8 mm        | ~ 60 mm       | 230 cm         | 1 piece |


Scale 1:1




#### With push handle, with isolated tip (Olympus HF connector)

| Article number | Working channel | Snare opening | Working length | PU      |
|----------------|-----------------|---------------|----------------|---------|
| 6522.401A-180  | Ø 2,8 mm        | ~ 40 mm       | 180 cm         | 1 piece |
| 6522.401A-230  | Ø 2,8 mm        | ~ 40 mm       | 230 cm         | 1 piece |
| 6522.601A-180  | Ø 2,8 mm        | ~ 60 mm       | 180 cm         | 1 piece |
| 6522.601A-230  | Ø 2,8 mm        | ~ 60 mm       | 230 cm         | 1 piece |

Scale 1:1





134°C

Autoclave

273°F


Reusable

# Hexagonal multifilament polypectomy snare



Without insulated tip

| Article number | Working channel | Snare opening | Working length | PU      |
|----------------|-----------------|---------------|----------------|---------|
| 37180/65-230   | Ø 2,0 mm        | - 25 mm       | 230 cm         | 1 piece |
| 37220/55-230   | Ø 2,8 mm        | - 20 mm       | 230 cm         | 1 piece |
| 37220/65-230   | Ø 2,8 mm        | - 25 mm       | 230 cm         | 1 piece |
| 37220/75-230   | Ø 2,8 mm        | - 35 mm       | 230 cm         | 1 piece |




Without insulated tip, with fixed handle (Olympus HF connector)

| Article number | Working channel | Snare opening | Working length | PU      |
|----------------|-----------------|---------------|----------------|---------|
| K37180A/65-230 | Ø 2,0 mm        | - 25 mm       | 230 cm         | 1 piece |
| K37220A/55-230 | Ø 2,8 mm        | - 20 mm       | 230 cm         | 1 piece |
| K37220A/65-230 | Ø 2,8 mm        | - 25 mm       | 230 cm         | 1 piece |
| K37220A/75-230 | Ø 2,8 mm        | - 35 mm       | 230 cm         | 1 piece |




Rotatable, without isolated tip

| Article number | Working channel | Snare opening | Working length | PU      |
|----------------|-----------------|---------------|----------------|---------|
| D37180/65-230  | Ø 2,0 mm        | - 25 mm       | 230 cm         | 1 piece |
| D37220/55-230  | Ø 2,8 mm        | - 20 mm       | 230 cm         | 1 piece |
| D37220/65-230  | Ø 2,8 mm        | - 25 mm       | 230 cm         | 1 piece |
| D37220/75-230  | Ø 2,8 mm        | - 35 mm       | 230 cm         | 1 piece |







| With insulated tip |                 |               |                |         |
|--------------------|-----------------|---------------|----------------|---------|
| Article number     | Working channel | Snare opening | Working length | PU      |
| 38180/65-230       | Ø 2,0 mm        | ~ 25 mm       | 230 cm         | 1 piece |
| 38220/55-230       | Ø 2,8 mm        | ~ 20 mm       | 230 cm         | 1 piece |
| 38220/65-230       | Ø 2,8 mm        | ~ 25 mm       | 230 cm         | 1 piece |
| 38220/75-230       | Ø 2,8 mm        | ~ 35 mm       | 230 cm         | 1 piece |



| With insulated tip, with fixed handle (Olympus HF connector) |                 |               |                |         |
|--|-----------------|---------------|----------------|---------|
| Article number   | Working channel | Snare opening | Working length | PU      |
| K38180A/65-230   | Ø 2,0 mm        | ~ 25 mm       | 230 cm         | 1 piece |
| K38220A/55-230   | Ø 2,8 mm        | ~ 20 mm       | 230 cm         | 1 piece |
| K38220A/65-230   | Ø 2,8 mm        | ~ 25 mm       | 230 cm         | 1 piece |
| K38220A/75-230   | Ø 2,8 mm        | ~ 35 mm       | 230 cm         | 1 piece |



| Rotatable, with insulated tip |                 |               |                |         |
|-------------------------------|-----------------|---------------|----------------|---------|
| Article number                | Working channel | Snare opening | Working length | PU      |
| D38180/65-230                 | Ø 2,0 mm        | ~ 25 mm       | 230 cm         | 1 piece |
| D38220/55-230                 | Ø 2,8 mm        | ~ 20 mm       | 230 cm         | 1 piece |
| D38220/65-230                 | Ø 2,8 mm        | ~ 25 mm       | 230 cm         | 1 piece |
| D38220/75-230                 | Ø 2,8 mm        | ~ 35 mm       | 230 cm         | 1 piece |



Reusable


Tube attachment for non-rotatable models



Reusable

Tube attachment for rotatable models

Reusable Handles




Reusable

## Handle with Olympus HF connector for tube attachment

Suitable for standard handle & handle for rotatable versions.


| Article number | Connection                          | PU      |
|----------------|-------------------------------------|---------|
| 6203LL         | Olympus HF connector for 4 mm cable | 1 piece |



Reusable

Fixed handle

# More products from Maslanka



Our product portfolio includes products for gastroenterology, cardiology and urology, as well as specialist solutions for surgery.

Our product catalogue contains a detailed breakdown of our products. Please do not hesitate to contact us if you require any further information.

Gastroenterology  
Cardiology  
Urology  
Neurology  
Customer-specific developments

# Reusable Polypectomy & mucosectomy snares



**H. + H. Maslanka**  
**Chirurgische Instrumente GmbH**  
Stockacherstraße 172  
78532 Tuttlingen / Germany

Tel.: +49 7461 9607 - 0  
Fax: +49 7461 - 47 32  
[info@maslanka.de](mailto:info@maslanka.de)  
[www.maslanka.de](http://www.maslanka.de)

CE 0483 07/2021 | EN