




Precise removal

Reusable polypectomy & mucosectomy snares

MASLANKA reusable polypectomy & mucosectomy snares

Benefit from the versatility of our snares




Symbol	Meaning
	Without insulated tip
	With insulated tip
	Multifilament
	Monofilament
	Rotatable

Flexibility & functionality in focus.

Our reusable snares, offer highest user comfort, whether with **multifilament or monofilament cutting wire**. The **non-insulated tip** enables selective coagulation of the removal site by the user after removal of polyps, without having to change instruments.

You can also benefit from our reusable **multifilament** snares **with insulated tip**. This prevents perforation during extension of the snare or coagulation during ablation of the polyp under poor visibility.




Reusable Symmetrical multifilament polypectomy snare

Our product for multiple use.



Without insulated tip

Article number	Working channel	Snare opening	Working length	PU
34181/55-230	Ø 2,0 mm	20 mm	230 cm	1 piece
34181/65-230	Ø 2,0 mm	25 mm	230 cm	1 piece
34181/75-230	Ø 2,0 mm	35 mm	230 cm	1 piece
34221/55-230	Ø 2,8 mm	20 mm	230 cm	1 piece
34221/65-230	Ø 2,8 mm	25 mm	230 cm	1 piece
34221/75-230	Ø 2,8 mm	35 mm	230 cm	1 piece




Without insulated tip, with fixed handle

Article number	Working channel	Snare opening	Working length	PU
K34221/55A-230	Ø 2,8 mm	20 mm	230 cm	1 piece
K34221/65A-230	Ø 2,8 mm	25 mm	230 cm	1 piece
K34221/75A-230	Ø 2,8 mm	35 mm	230 cm	1 piece




Rotatable, without insulated tip


Article number	Working channel	Snare opening	Working length	PU
D34221/55-230	Ø 2,8 mm	20 mm	230 cm	1 piece
D34221/65-230	Ø 2,8 mm	25 mm	230 cm	1 piece
D34221/75-230	Ø 2,8 mm	35 mm	230 cm	1 piece





With insulated tip				
Article number	Working channel	Snare opening	Working length	PU
34180/55-230	Ø 2,0 mm	20 mm	230 cm	1 piece
34180/65-230	Ø 2,0 mm	25 mm	230 cm	1 piece
34180/75-230	Ø 2,0 mm	35 mm	230 cm	1 piece
34220/55-230	Ø 2,8 mm	20 mm	230 cm	1 piece
34220/65-230	Ø 2,8 mm	25 mm	230 cm	1 piece
34220/75-230	Ø 2,8 mm	35 mm	230 cm	1 piece



With insulated tip, with fixed handle (Olympus HF connector)				
Article number	Working channel	Snare opening	Working length	PU
K34220/55A-230	Ø 2,8 mm	20 mm	230 cm	1 piece
K34220/65A-230	Ø 2,8 mm	25 mm	230 cm	1 piece
K34220/75A-230	Ø 2,8 mm	35 mm	230 cm	1 piece



Rotatable, with insulated tip				
Article number	Working channel	Snare opening	Working length	PU
D34220/55-230	Ø 2,8 mm	20 mm	230 cm	1 piece
D34220/65-230	Ø 2,8 mm	25 mm	230 cm	1 piece
D34220/75-230	Ø 2,8 mm	35 mm	230 cm	1 piece




Reusable



Symmetrical monofilament polypectomy snare

Our product for multiple use.



Without insulated tip				
Article number	Working channel	Snare opening	Working length	PU
34221M/35-230	Ø 2,8 mm	10 mm	230 cm	1 piece
34221M/55-230	Ø 2,8 mm	20 mm	230 cm	1 piece
34221M/65-230	Ø 2,8 mm	25 mm	230 cm	1 piece
34221M/75-230	Ø 2,8 mm	30 mm	230 cm	1 piece





Reusable

Symmetrical multifilament mucosectomy snare


Our product for multiple use. Available only with fixed handle.



Without insulated tip (Olympus HF connector)

Article number	Working channel	Snare opening	Working length	PU
K34226A/75-230	Ø 2,8 mm	30 mm	230 cm	1 piece

Scale 1:1



Reusable

Asymmetrical multifilament polypectomy snare

Our product for multiple use.

Without insulated tip				
Article number	Working channel	Snare opening	Working length	PU
3218-180	Ø 2,0 mm	15 mm	160 cm	1 piece
3218-230	Ø 2,0 mm	15 mm	230 cm	1 piece
3322-180	Ø 2,8 mm	25 mm	160 cm	1 piece
3322-230	Ø 2,8 mm	25 mm	230 cm	1 piece

Scale 1:1

With insulated tip				
Article number	Working channel	Snare opening	Working length	PU
3022-230	Ø 2,8 mm	15 mm	230 cm	1 piece
3122-230	Ø 2,8 mm	25 mm	230 cm	1 piece

Scale 1:1

Without insulated tip, with fixed handle (Olympus HF connector)				
Article number	Working channel	Snare opening	Working length	PU
K3218A-180	Ø 2,0 mm	15 mm	160 cm	1 piece
K3218A-230	Ø 2,0 mm	15 mm	230 cm	1 piece
K3322A-180	Ø 2,8 mm	25 mm	160 cm	1 piece
K3322A-230	Ø 2,8 mm	25 mm	230 cm	1 piece


With insulated tip, with fixed handle (Olympus HF connector)				
Article number	Working channel	Snare opening	Working length	PU
K3022A-230	Ø 2,8 mm	15 mm	230 cm	1 piece
K3122A-230	Ø 2,8 mm	25 mm	230 cm	1 piece

Rotatable, without isolated tip				
Article number	Working channel	Snare opening	Working length	PU
D3322-230	Ø 2,8 mm	25 mm	230 cm	1 piece

Scale 1:1

Rotatable with insulated tip				
Article number	Working channel	Snare opening	Working length	PU
D3122-230	Ø 2,8 mm	25 mm	230 cm	1 piece

Scale 1:1



Reusable



Asymmetrical multifilament polypectomy snare

Our product for multiple use.




With push handle, without isolated tip (Olympus HF connector)				
Article number	Working channel	Snare opening	Working length	PU
6522.40A-180	Ø 2,8 mm	40 mm	180 cm	1 piece
6522.40A-230	Ø 2,8 mm	40 mm	230 cm	1 piece
6522.60A-180	Ø 2,8 mm	60 mm	180 cm	1 piece
6522.60A-230	Ø 2,8 mm	60 mm	230 cm	1 piece


Scale 1:1



With push handle, with isolated tip (Olympus HF connector)				
Article number	Working channel	Snare opening	Working length	PU
6522.401A-180	Ø 2,8 mm	40 mm	180 cm	1 piece
6522.401A-230	Ø 2,8 mm	40 mm	230 cm	1 piece
6522.601A-180	Ø 2,8 mm	60 mm	180 cm	1 piece
6522.601A-230	Ø 2,8 mm	60 mm	230 cm	1 piece

Scale 1:1





Reusable


Hexagonal multifilament polypectomy snare

Our product for multiple use.



Without insulated tip

Article number	Working channel	Snare opening	Working length	PU
37180/65-230	Ø 2,0 mm	25 mm	230 cm	1 piece
37220/55-230	Ø 2,8 mm	20 mm	230 cm	1 piece
37220/65-230	Ø 2,8 mm	25 mm	230 cm	1 piece
37220/75-230	Ø 2,8 mm	35 mm	230 cm	1 piece




Without insulated tip, with fixed handle (Olympus HF connector)

Article number	Working channel	Snare opening	Working length	PU
K37180A/65-230	Ø 2,0 mm	25 mm	230 cm	1 piece
K37220A/55-230	Ø 2,8 mm	20 mm	230 cm	1 piece
K37220A/65-230	Ø 2,8 mm	25 mm	230 cm	1 piece
K37220A/75-230	Ø 2,8 mm	35 mm	230 cm	1 piece



Rotatable, without isolated tip


Article number	Working channel	Snare opening	Working length	PU
D37180/65-230	Ø 2,0 mm	25 mm	230 cm	1 piece
D37220/55-230	Ø 2,8 mm	20 mm	230 cm	1 piece
D37220/65-230	Ø 2,8 mm	25 mm	230 cm	1 piece
D37220/75-230	Ø 2,8 mm	35 mm	230 cm	1 piece




With insulated tip				
Article number	Working channel	Snare opening	Working length	PU
38180/65-230	Ø 2,0 mm	25 mm	230 cm	1 piece
38220/55-230	Ø 2,8 mm	20 mm	230 cm	1 piece
38220/65-230	Ø 2,8 mm	25 mm	230 cm	1 piece
38220/75-230	Ø 2,8 mm	35 mm	230 cm	1 piece



With insulated tip, with fixed handle (Olympus HF connector)				
Article number	Working channel	Snare opening	Working length	PU
K38180A/65-230	Ø 2,0 mm	25 mm	230 cm	1 piece
K38220A/55-230	Ø 2,8 mm	20 mm	230 cm	1 piece
K38220A/65-230	Ø 2,8 mm	25 mm	230 cm	1 piece
K38220A/75-230	Ø 2,8 mm	35 mm	230 cm	1 piece




Rotatable, with insulated tip				
Article number	Working channel	Snare opening	Working length	PU
D38180/65-230	Ø 2,0 mm	25 mm	230 cm	1 piece
D38220/55-230	Ø 2,8 mm	20 mm	230 cm	1 piece
D38220/65-230	Ø 2,8 mm	25 mm	230 cm	1 piece
D38220/75-230	Ø 2,8 mm	35 mm	230 cm	1 piece



Reusable

Tube attachment for non-rotatable models



Reusable

Tube attachment for rotatable models

Our product for multiple use.

Reusable Handles




Reusable

Handle with Olympus HF connector for tube attachment

Our product for multiple use. Suitable for standard handle & handle for rotatable versions.


Article number	Connection	PU
6203LL	Olympus HF connector for 4 mm cable	1 piece



Reusable

Fixed handle

More products from Maslanka



Our product portfolio includes products for gastroenterology, cardiology and urology, as well as specialist solutions for surgery.

Our product catalogue contains a detailed breakdown of our products. Please do not hesitate to contact us if you require any further information.

Gastroenterology
Cardiology
Urology
Neurology
Customer-specific developments

Reusable Polypectomy & mucosectomy snares



H. + H. Maslanka
Chirurgische Instrumente GmbH
Stockacherstraße 172
78532 Tuttlingen / Germany

Tel.: +49 7461 9607 - 0
Fax: +49 7461 - 47 32
info@maslanka.de
www.maslanka.de

05/2020 | EN

CE 0483