



Securely held

 Cystoscopy

MASLANKA Reusable urology foreign body forceps


Benefit from the advantages of our forceps

Luer Lock connector for cleaning the inner lumen


Force limitation

From the handle to the coordinated angle – the entire product is ergonomic and therefore offers comfortable use

For finger support



Heat shrink tubing increases the kink protection stability at the proximal end




Rigid = tube



Semi-flexible = Thin-walled tube



Flexible = Spiral and monofilament pull wire





Reusable

Interlocking foreign body forceps, 4.6 mm long

Flexible: spiral ground over the entire length

Article number	Instrument diameter	Opening width	Volumes	Mouthpiece length	Working length	PU
S04616-040	1,6 mm	- 8,4 mm	- 5,2 mm ³	- 4,6 mm	40 cm	1 x
S04616-060	1,6 mm	- 8,4 mm	- 5,2 mm ³	- 4,6 mm	60 cm	1 x

Rigid: tube and wire as inner cable

Article number	Instrument diameter	Opening width	Volumes	Mouthpiece length	Working length	PU
RS04616-060	1,6 mm	- 8,4 mm	- 5,2 mm ³	- 4,6 mm	60 cm	1 x




Reusable

Foreign body forceps interlocking teeth, 6.9 mm long

Rigid: tube and wire as inner cable

Article number	Instrument diameter	Opening width	Volumes	Mouthpiece length	Working length	PU
RS04316-040	1,6 mm	- 11,2 mm	- 3 mm ³	- 6,9 mm	40 cm	1 x
RS04316-060	1,6 mm	- 11,2 mm	- 3 mm ³	- 6,9 mm	60 cm	1 x



Reusable

Interlocking foreign body forceps

Flexible: spiral ground over the entire length

Article number	Instrument diameter	Opening width	Mouthpiece length	Working length	PU
S05212-040	1,2 mm	- 6 mm	- 3,8 mm ³	40 cm	1 x




Reusable

Interlocking foreign body grasping forceps, one-sided

Flexible: spiral ground over the entire length

Article number	Instrument diameter	Opening width	Volumes	Mouthpiece length	Working length	PU
S05316-060	1,6 mm	- 5,2 mm	- 5,3 mm ³	- 3,7 mm	60 cm	1 x



Reusable

Foreign body grasping forceps, alligator jaws

Flexible: spiral ground over the entire length

Article number	Instrument diameter	Opening width	Mouthpiece length	Working length	PU
S1816-040	1,6 mm	- 5,5 mm	- 3,7 mm	40 cm	1 x

Rigid: tube and wire as inner cable

Article number	Instrument diameter	Opening width	Mouthpiece length	Working length	PU
RS1816-060	1,6 mm	- 5,5 mm	- 3,7 mm	60 cm	1 x
RS1822-040	2,3 mm	- 8,5 mm	- 5,7 mm	40 cm	1 x



Reusable

Foreign body forceps for steinstrasse

Flexible: spiral ground over the entire length

Article number	Instrument diameter	Opening width	Volumes	Mouthpiece length	Working length	PU
S04116-030	1,6 mm	- 13,8 mm	- 10,5 mm ³	- 8,2 mm	30 cm	1 x

Rigid: tube and wire as inner cable

Article number	Instrument diameter	Opening width	Volumes	Mouthpiece length	Working length	PU
RS04116-060	1,6 mm	- 13,8 mm	- 10,5 mm ³	- 8,2 mm	60 cm	1 x



Reusable DJ forceps

Flexible: spiral ground over the entire length

Article number	Instrument diameter	Opening width	Mouthpiece length	Working length	PU
S08322-040	2,3 mm	- 8,0 mm	- 6,1 mm	40 cm	1 x
S08322-060	2,3 mm	- 8,0 mm	- 6,1 mm	60 cm	1 x

Rigid: tube and wire as inner cable

Article number	Instrument diameter	Opening width	Mouthpiece length	Working length	PU
RS08322-040	2,3 mm	- 8,0 mm	- 6,1 mm	40 cm	1 x
RS08322-060	2,3 mm	- 8,0 mm	- 6,1 mm	60 cm	1 x





Reusable

Rat-tooth grasping forceps

Flexible: spiral ground over the entire length

Article number	Instrument diameter	Opening width	Mouthpiece length	Working length	PU
S08416-040	1,6 mm	- 9,0 mm	- 4,2 mm	40 cm	1 x
S08416-060	1,6 mm	- 9,0 mm	- 4,2 mm	60 cm	1 x
S08422-040	2,3 mm	- 5,6 mm	- 4,2 mm	40 cm	1 x
S08422-060	2,3 mm	- 5,6 mm	- 4,2 mm	60 cm	1 x

More products from Maslanka



Products for gastroenterology, cardiology and urology as well as special solutions for surgery are part of our product range.

You will find a detailed overview of our products in our product catalogue. Please do not hesitate to contact us if you have any further questions.

Gastroenterology
Cardiology
Urology
Neurology
Customised developments

Cystoscopy

Reusable urology foreign body forceps



H. + H. Maslanka
Chirurgische Instrumente GmbH
Stockacherstraße 172
78532 Tuttlingen / Germany

Tel.: +49 (0) 7461 9607 - 0
Fax: +49 (0) 7461 - 47 32
info@maslanka.de
www.maslanka.de

03/2022 | EN

CE 0483