



Disposable





# Precise removal

Disposable polypectomy & mucosectomy snares

MASLANKA disposable polypectomy & mucosectomy snares

# Benefit from the advantages of our snares



## Flexibility & functionality in focus.

Our **multifilament** slings offer highest user comfort for single use. The **non-insulated tip** enables selective coagulation of the removal site by the user after removal of polyps, without changing instruments.


 Disposable

## Symmetrical multifilament polypectomy snare

Our product for single use. Available only with fixed-mounted handle.



Article number	Working channel	Snare opening	Working length	PU
E34221A/55-180	Ø 2,8 mm	20 mm	180 cm	10 piece
E34221A/55-230	Ø 2,8 mm	20 mm	230 cm	10 piece
E34221A/65-180	Ø 2,8 mm	25 mm	180 cm	10 piece
E34221A/65-230	Ø 2,8 mm	25 mm	230 cm	10 piece
E34221A/75-180	Ø 2,8 mm	30 mm	180 cm	10 piece
E34221A/75-230	Ø 2,8 mm	30 mm	230 cm	10 piece




 Disposable

## Symmetrical multifilament mucosectomy snare

Our product for single use. Available only with fixed-mounted handle.



Article number	Working channel	Snare opening	Working length	PU
E34224A/65-180	Ø 2,8 mm	25 mm	180 cm	10 piece
E34224A/65-230	Ø 2,8 mm	25 mm	230 cm	10 piece
E34224A/75-180	Ø 2,8 mm	30 mm	180 cm	10 piece
E34224A/75-230	Ø 2,8 mm	30 mm	230 cm	10 piece




 Disposable

## Symmetrical monofilament polypectomy snare

Our product for single use. Available only with fixed-mounted handle.



Article number	Working channel	Snare opening	Working length	PU
E34221MA/35-180	Ø 2,8 mm	10 mm	180 cm	10 piece
E34221MA/35-230	Ø 2,8 mm	10 mm	230 cm	10 piece
E34221MA/55-180	Ø 2,8 mm	20 mm	180 cm	10 piece
E34221MA/55-230	Ø 2,8 mm	20 mm	230 cm	10 piece
E34221MA/75-180	Ø 2,8 mm	30 mm	180 cm	10 piece
E34221MA/75-230	Ø 2,8 mm	30 mm	230 cm	10 piece




 Disposable

## Hexagonal multifilament polypectomy snare

Our product for single use. Available only with fixed-mounted handle.



Article number	Working channel	Snare opening	Working length	PU
E34221A/55-180	Ø 2,8 mm	15 mm	180 cm	10 piece
E34221A/55-230	Ø 2,8 mm	15 mm	230 cm	10 piece
E34221A/65-180	Ø 2,8 mm	25 mm	180 cm	10 piece
E34221A/65-230	Ø 2,8 mm	25 mm	230 cm	10 piece
E34221A/75-180	Ø 2,8 mm	30 mm	180 cm	10 piece
E34221A/75-230	Ø 2,8 mm	30 mm	230 cm	10 piece




 Disposable

## Asymmetrical multifilament polypectomy snare

Our product for single use. Available only with fixed-mounted handle.




Article number	Working channel	Snare opening	Working length	PU
E3222A-180	Ø 2,8 mm	15 mm	180 cm	10 piece
E3222A-230	Ø 2,8 mm	15 mm	230 cm	10 piece
E3322A-180	Ø 2,8 mm	25 mm	180 cm	10 piece
E3322A-230	Ø 2,8 mm	25 mm	230 cm	10 piece



## Serrated snare

with sharp teeth for better tissue grip




## Disposable

### Symmetrical serrated mucosectomy snare

Our product for single use. Available only with fixed-mounted handle.




Article number	Working channel	Snare opening	Working length	PU
E342213A/55-180	Ø 2,8 mm	20 mm	180 cm	10 piece
E342213A/55-230	Ø 2,8 mm	20 mm	230 cm	10 piece
E342213A/65-180	Ø 2,8 mm	25 mm	180 cm	10 piece
E342213A/65-230	Ø 2,8 mm	25 mm	230 cm	10 piece
E342213A/75-180	Ø 2,8 mm	30 mm	180 cm	10 piece
E342213A/75-230	Ø 2,8 mm	30 mm	230 cm	10 piece



## With thermal coating

for better conductivity of the RF current and a smoother cutting surface on the tissue through very low coagulation.




## Disposable


### Symmetrical multifilament polypectomy snare with thermo-layer

Our product for single use. Available only with fixed-mounted handle.



Article number	Working channel	Snare opening	Working length	PU
E34223A/55-230	Ø 2,8 mm	20 mm	230 cm	10 piece
E34223A/65-230	Ø 2,8 mm	25 mm	230 cm	10 piece
E34223A/75-230	Ø 2,8 mm	30 mm	230 cm	10 piece





Ergonomic handle  
for **better handling**

RF connection

Due to its fixed handle, the snare  
is **ready for use immediately** after  
plugging in the RF cable

## **Disposable** Fixed handle

Our product for single use.

# Disposable Polypectomy & mucosectomy snares



**H. + H. Maslanka**  
**Chirurgische Instrumente GmbH**  
Stockacherstraße 172  
78532 Tuttlingen / Germany

Tel.: +49 7461 9607 - 0  
Fax: +49 7461 - 47 32  
info@maslanka.de  
www.maslanka.de

04/2020 | EN

**CE 0483**