



Disposable





Precise removal

Disposable polypectomy & mucosectomy snares

MASLANKA disposable polypectomy & mucosectomy snares

Benefit from the advantages of our snares



Flexibility & functionality in focus.

Our **multifilament** slings offer highest user comfort for single use.
The **non-insulated tip** enables selective coagulation of the removal site by the user after removal of polyps, without changing instruments.


ⓧ Disposable

Symmetrical multifilament polypectomy snare

Available only with fixed-mounted handle.



Article number	Working channel	Snare opening	Working length	PU
E34221A/35-180	Ø 2,8 mm	~ 10 mm	180 cm	10 pieces
E34221A/35-230	Ø 2,8 mm	~ 10 mm	230 cm	10 pieces
E34221A/45-180	Ø 2,8 mm	~ 15 mm	180 cm	10 pieces
E34221A/45-230	Ø 2,8 mm	~ 15 mm	230 cm	10 pieces
E34221A/55-180	Ø 2,8 mm	~ 20 mm	180 cm	10 pieces
E34221A/55-230	Ø 2,8 mm	~ 20 mm	230 cm	10 pieces
E34221A/65-180	Ø 2,8 mm	~ 25 mm	180 cm	10 pieces
E34221A/65-230	Ø 2,8 mm	~ 25 mm	230 cm	10 pieces
E34221A/75-180	Ø 2,8 mm	~ 35 mm	180 cm	10 pieces
E34221A/75-230	Ø 2,8 mm	~ 35 mm	230 cm	10 pieces




ⓧ Disposable

Symmetrical monofilament polypectomy snare

Available only with fixed-mounted handle.



Article number	Working channel	Snare opening	Working length	PU
E34221MA/35-180	Ø 2,8 mm	~ 10 mm	180 cm	10 pieces
E34221MA/35-230	Ø 2,8 mm	~ 10 mm	230 cm	10 pieces
E34221MA/45-180	Ø 2,8 mm	~ 15 mm	180 cm	10 pieces
E34221MA/45-230	Ø 2,8 mm	~ 15 mm	230 cm	10 pieces
E34221MA/55-180	Ø 2,8 mm	~ 20 mm	180 cm	10 pieces
E34221MA/55-230	Ø 2,8 mm	~ 20 mm	230 cm	10 pieces
E34221MA/75-180	Ø 2,8 mm	~ 30 mm	180 cm	10 pieces
E34221MA/75-230	Ø 2,8 mm	~ 30 mm	230 cm	10 pieces




ⓧ Disposable

Asymmetrical multifilament polypectomy snare

Available only with fixed-mounted handle.



Article number	Working channel	Snare opening	Working length	PU
E3222A-180	Ø 2,8 mm	~ 15 mm	180 cm	10 pieces
E3222A-230	Ø 2,8 mm	~ 15 mm	230 cm	10 pieces
E3322A-180	Ø 2,8 mm	~ 25 mm	180 cm	10 pieces
E3322A-230	Ø 2,8 mm	~ 25 mm	230 cm	10 pieces




ⓧ Disposable

Symmetrical multifilament mucosectomy snare

Available only with fixed-mounted handle.



Article number	Working channel	Snare opening	Working length	PU
E34224A/65-180	Ø 2,8 mm	~ 25 mm	180 cm	10 pieces
E34224A/65-230	Ø 2,8 mm	~ 25 mm	230 cm	10 pieces
E34224A/75-180	Ø 2,8 mm	~ 30 mm	180 cm	10 pieces
E34224A/75-230	Ø 2,8 mm	~ 30 mm	230 cm	10 pieces




ⓧ Disposable

Hexagonal multifilament polypectomy snare

Available only with fixed-mounted handle.




Article number	Working channel	Snare opening	Working length	PU
E37221A/55-180	Ø 2,8 mm	~ 20 mm	180 cm	10 pieces
E37221A/55-230	Ø 2,8 mm	~ 20 mm	230 cm	10 pieces
E37221A/65-180	Ø 2,8 mm	~ 25 mm	180 cm	10 pieces
E37221A/65-230	Ø 2,8 mm	~ 25 mm	230 cm	10 pieces
E37221A/75-180	Ø 2,8 mm	~ 35 mm	180 cm	10 pieces
E37221A/75-230	Ø 2,8 mm	~ 35 mm	230 cm	10 pieces



Serrated snare

with sharp teeth for better tissue grip




 Disposable

Symmetrical serrated mucosectomy snare

Available only with fixed-mounted handle.




Article number	Working channel	Snare opening	Working length	PU
E342213A/55-180	Ø 2,8 mm	- 20 mm	180 cm	10 pieces
E342213A/55-230	Ø 2,8 mm	- 20 mm	230 cm	10 pieces
E342213A/65-180	Ø 2,8 mm	- 25 mm	180 cm	10 pieces
E342213A/65-230	Ø 2,8 mm	- 25 mm	230 cm	10 pieces
E342213A/75-180	Ø 2,8 mm	- 35 mm	180 cm	10 pieces
E342213A/75-230	Ø 2,8 mm	- 35 mm	230 cm	10 pieces



With thermal coating

for better conductivity of the RF current and a smoother cutting surface on the tissue through very low coagulation.




Disposable


Symmetrical multifilament polypectomy snare with thermo-layer

Available only with fixed-mounted handle.



Article number	Working channel	Snare opening	Working length	PU
E34223A/55-230	Ø 2,8 mm	- 20 mm	230 cm	10 pieces
E34223A/65-230	Ø 2,8 mm	- 25 mm	230 cm	10 pieces
E34223A/75-230	Ø 2,8 mm	- 35 mm	230 cm	10 pieces





Disposable Fixed handle

Disposable Polypectomy & mucosectomy snares



H. + H. Maslanka
Chirurgische Instrumente GmbH
Stockacherstraße 172
78532 Tuttlingen / Germany

Tel.: +49 7461 9607 - 0
Fax: +49 7461 - 47 32
info@maslanka.de
www.maslanka.de

06/2021 | EN

CE 0483